

Uniaxial Force Plate

Brandon Segal, Bryanna Wands, Lubaina Shakir, Mayank Parmar




College of Engineering

Spinal Cord Injury

Annual incidence - 54 cases/million

Approximately 300,000 persons currently

Less than 1 percent experience complete recovery



Current Research

Electrical Stimulation


Optogenetics




Magnetogenetics



Chemogenetics - DREADDs





Proposed Solution – Force Plate

- Tool for kinetic analysis of rodent gait
- Appropriate for rodent model
- Cost-efficient
- Avoids inaccuracies due to loose skin



Design Requirements of force plate

Objective

- Have a linear response in the range of forces
- Reasonable operable procedure

Numerical

- Have a natural frequency $> 100\text{Hz}$
- Sensitive to weights as low as 0.1 N
- Carry loads up to 5 N

Technical

- Have a response independent of where the force is exerted
- Have a sufficient coefficient of safety to protect the specimen

Animal Care

Mayank Parmar

Training Inputs

Handling

- Must allow Rats to acclimate to human interaction

Environment

- Rats must acclimate to treadmill environment

Responsiveness

- Rats must be trained to respond to food reward system



Acquired Outputs


Cooperation

- Rats acclimated to protocol procedures requiring gentle handling.
- Rats have become actively comfortable in the hands of the caretaker and handler.

Performance

- Rats adjusted to being placed in the treadmill enclosure.
- At various speeds, the rats were able to maintain their speed and avoid major collisions with the enclosure.

Response to Reward System

- Rats were trained to receive food rewards after runs in order to give them incentive to perform the action.
 - Food reward system proved to display the calmness of the rats on the treadmill between speed intervals
- 

Example Training Schedule

Training Speed(cm/sec)	Run Time (sec)	Rest Time Between Intervals(sec)
16	30-60	90-120
20	30-60	90-120
24	30-60	90-120
28	30-60	90-120
32	30-60	90-120

WEEK 1

Day 1



Example Training Schedule (continued)

Training Speed(cm/sec)	Run Time (sec)	Rest Time Between Intervals(sec)
16	60-120	60-90
20	60-120	60-90
24	60-120	60-90
28	60-120	60-90
32	60-120	60-90

WEEK 1

Day 2/3



Example Training Schedule (continued)

Training Speed(cm/sec)	Run Time (sec)	Rest Time Between Intervals(sec)
16	300	60
20	300	60
24	300	60
28	300	60-
32	300	60

WEEK 2-6

Day 1-3




Questionable Aspects of Training Protocol

- Why training occurs?
- Why these specific speeds are chosen?
- Why these specific time intervals?



Mechanical Requirements for the force plate



Sensors

Bryanna Wands

Design Inputs

Sensitivity

- Must sense forces 0-5N(Zumwalt)

Linearity

- Data must show linear force versus strain relationship

Stability


- Must be stable over time and changing conditions



Design Outputs

Strain Gauges

- Simple circuitry
- Half-bridge configuration
 - Increases sensitivity
 - Negates effects of temperature change



Evolution of Sensor Design

Force

- Piezoelectric: Voltage proportional to applied force using piezoelectric effect
- Loadcells: Manufactured sensor packages measuring compression forces

Displacement

- Hall Effect Sensors: Assembly of springs and magnets measuring the change in magnetic field due to the magnet's deflection
- Optical Sensors: Use Phototransistors to measure change in lumens due to the deflection of a LED or Laser

Strain

- Strain Gauges: Measure strain using the deformation of metal elements in a sensor
- 

Evolution of Sensor Design

Sensors		Pros	Cons
Piezoelectric	1	<ul style="list-style-type: none">• Directly Measures Force• Simple Circuitry	<ul style="list-style-type: none">• Can Only measure force in Z direction• Can only measure dynamic forces
Load Cell		<ul style="list-style-type: none">• Directly Measures Force	<ul style="list-style-type: none">• Prohibitively expensive• Only measures force in Z direction
HES		<ul style="list-style-type: none">• Displacement is linear to force• No need for opamp	<ul style="list-style-type: none">• Difficult to manufacture• Difficult to avoid crosstalk
Phototransistors	2	<ul style="list-style-type: none">• Not Temperature Dependent• Had Existing Drivers Available	<ul style="list-style-type: none">• Difficult to Align• Impossible to implement elsewhere
Strain Gauges	3	<ul style="list-style-type: none">• Proven Design• Able to Implement elsewhere	<ul style="list-style-type: none">• Temperature Dependent• Strain Gauges may detach.

Evolution of Sensor Design


Tonia Hsieh Phototransistor
Driver:

Inputs(V): CH1INA, CH1INB

Outputs(V): CH1OUTA,
CH1OUTB, CH1OUTAB


Components:

- 2 Low Pass filters amplifiers
- 15 V in
- LED Power
- BNC outputs to DAQ



Evolution of Sensor Design

- Strain Gauge Circuit
- Half-Bridge Configuration
- Measure the Voltage over R_g






Geometry & Materials

Brandon Segal

Geometry Design Inputs



- SolidWorks
- Rat Sizes
- Heglund Cantilever Model
- Hseih Box Model

Evolution of Geometry

Platform

- This was theorized when developing ideas on how to incorporate piezoelectric sensors

Cantilever


- Adaption to the Heglund model created in 1981
- Allows for multiple dimensions of force

Evolution of Geometry

The box acts as two
spring blades

A fusion of Tonia
Hsieh's and
Heglund's Designs.

Measure strain or
displacement.



Geometry Design Outputs

Allows for
sensor
readings

Proven Model

Testable
Dimensions

Key Values can
be found using
basic equations

Allows Future
Work



Material Design Inputs




Strain Values

Displacement

Natural Frequency

Stress Values

- Agreed upon Force plate Geometry
 - Equations to calculate key values
 - Deflection.jl Code
 - SolidWorks Model
- 

Evolution of Material Choice

Copper

- Pliable and machinable metal used in highly sensitive force plates (Hsieh 2006)

ABS Plastic



- Acrylonitrile butadiene styrene- Durable plastic available for 3D printing at Temple

PLC Plastic

- Polycaprolactone – Fragile plastic used by most 3D printers

Evolution of Material Choice

EPSX



Youngs Modulus
for ABS Plastic: 2.00 GPa

Average Strain Values at Box
edge (tension side): 37.5 μ


Calculated Strain Values: 36.
75 μ

Material Design Outputs

Spring Blade Dimensions


- Able to finalize Solidworks dimensions
- Output results to CSV file
- Analyzed Results to pick the final dimensions

SolidWorks Simulation

- Used to validate hand calculations
 - Average Values correlated within 5% of expected values
 - Dimensions could hold chosen sensors
- 

Testing of the Equipment

Follow A Similar procedure as Zumwalt



Force Plate Operation

Lubaina Shakir

Experimental Protocol


Normal Baseline

- Healthy rats


SCI Baseline

- Forces are captured after hemisection of spinal cord


DREADDs + CNO

- Rats transfected with DREADDs constructed as per protocol
 - CNO given just prior to data collection
- 

Model for Overground Course




Data Collection Output (Theoretical)



Zumwalt (2005)

Ensuring Performance of the Equipment




Ensuring Animal Safety




Animal Training

- Make sure all personnel are trained to handle animals



Equipment Readiness

- Ensure that Surfaces are flush with the course



Structural Integrity

- Implement a Factor of safety when choosing force plate geometries



Sustainability and Future Use

- Extra Strain Gauges will be donated to the Spence Lab
- Life of Strain Gauge exceeds 10,000,000 cycles at small loads
- The structure will remain in the Spence Lab for future teams to use as a reference
- Designs and documents will be uploaded to Github
- Make the design open source solution for small labs





Questions?